


Reusable Suction Cups



Ref. 077.0177 (Ø 40 mm)
 Ref. 077.0176 (Ø 50 mm)
 Ref. 077.0175 (Ø 60 mm)

Malmström Cup

- Autoclavable stainless steel
- Used when the fetus is in the occiputo-anterior position
- Cup design allows for good seal to help minimise pop-offs
- Sizes: Ø 40, 50 and 60 mm



Ref. 077.0173 (Ø 40 mm)
 Ref. 077.0172 (Ø 50 mm)
 Ref. 077.0171 (Ø 60 mm)

Bird Cup

- Autoclavable stainless steel
- Based on the design of the Malmström Cup, with separate suction and traction ports
- Sizes: Ø 40, 50 or 60 mm



Ref. 077.0174

Bird Posterior Cup

- Autoclavable stainless steel
- Especially designed for use when the fetus is in the occiputo-posterior position
- Vacuum port on side of suction cup and low profile to ensure easier insertion
- Size: Ø 50 mm



Ref. 077.0761 (Ø 50 mm) blue
 Ref. 077.0762 (Ø 60 mm) yellow

Silc Cup

- Soft suction cup for vacuum assisted delivery
- With silicone tubing, 1.5 m and 1 coupling piece
- Sizes: Ø 50 and 60 mm



Ref. 077.0078

Silc Cup with trumpet valve

- Attach gently to fetal scalp
- One-piece design for easy handling and use
- Smooth surface with longitudinal ridge makes it possible to observe rotation
- Size: Ø 60 mm



Ref. 077.0080

CaesAid Cup for Caesarean sections

- Autoclavable medical grade silicone
- Gentle support for delivering fetal head during Caesarean sections
- Can be attached to the pump used for suctioning during the operation.
- Size: Ø 60 mm


Disposable Suction Cups



Ref. 077.0791

Single Use Silc Cup

- Attaches gently to fetal scalp
- Smooth surface with longitudinal ridge makes it possible to observe rotation
- Sterile with 150 cm tubing
- Size: Ø 60 mm



Ref. 077.0792

Single Use Bird Cup

- Polycarbonate version of the reusable Bird Cup
- Traction handle can be tilted for use when the fetus is in the occiputo-anterior or posterior position
- Sterile with 150 cm tubing
- Size: Ø 50 mm